


I was able to successfully display the screen below using this URL !!!

Conquest DICOM server - version 1.5.0d

127.0.0.1/app/newweb/

Welcome to Conquest DICOM server - version 1.5.0d

Written by Conquest Project



DICOM server 'PacsPlus' (version 1.5.0d, port 4810, bits 64) was started on Wed Nov 01 22:29:23 2023
Old JPEG decoder=0, JPEGLIB jpeg codec=1, LIBJASPER jpeg2000 codec=1
Run time (s) total 2, query 0, load 0, save 0, compress 0, process 0, gpps 0
Associations=0; Threads=85 (1 open); Images sent=0, recieved=0, forwarded=0
Images printed=0, in color=0
Activity: Echo:0, Find:0, Move:0, Unknown:0, gpps:1310
Images (de)compressed: NKI 0, JPEG 0, JPEG2000 0, RLE 0, Planes removed 0, Palettes removed 0, Downsize 0
Space on MAGD : 17564 MByte
Database type: ODBC connection

Local Patient List Patient ID: Name:

Local Studies List Patient ID: Name: Date:

Upload file to enter into server (dcm/v2/HL7/zip/7z/gz/tar): ファイルが選択されていません


When I searched by patient ID and pressed the {submit} button, an error message was displayed.

Conquest DICOM server - version 1.5.0d

127.0.0.1/app/newweb/

Welcome to Conquest DICOM server - version 1.5.0d

Written by Conquest Project



DICOM server 'PacsPlus' (version 1.5.0d, port 4810, bits 64) was started on Wed Nov 01 22:29:23 2023
Old JPEG decoder=0, JPEGLIB jpeg codec=1, LIBJASPER jpeg2000 codec=1
Run time (s) total 2, query 0, load 0, save 0, compress 0, process 0, gpps 0
Associations=0; Threads=85 (1 open); Images sent=0, recieved=0, forwarded=0
Images printed=0, in color=0
Activity: Echo:0, Find:0, Move:0, Unknown:0, gpps:1310
Images (de)compressed: NKI 0, JPEG 0, JPEG2000 0, RLE 0, Planes removed 0, Palettes removed 0, Downsize 0
Space on MAGD : 17564 MByte
Database type: ODBC connection

Local Patient List Patient ID: Name:

Local Studies List Patient ID: Name: Date:

Upload file to enter into server (dcm/v2/HL7/zip/7z/gz/tar): ファイルが選択されていません

127.0.0.1/app/newweb/?mode=liststudies&key=&patientidmatch=420137&patien...

Welcome to Conquest DICOM server - version 1.5.0d

*** lua run error liststudies.lua:86: attempt to concatenate a nil value in 'dofile('liststudies.lua')'